Interpretation Guide



Pharmacogenomic (PGx) Report - for your healthcare provider

The following PGx report is a clinical decision support tool based on individual genetic results. It contributes to a better understanding and prediction of medication response and tolerability. This test does not predict the risk of any health problem. Since response to medications is multifactorial, clinical judgment supersedes any recommendations provided.

The report notifies you if the patient carries any genetic variant that can alter the following pharmacological parameters:

- pharmacokinetics: overall **exposure** to a medication depending on metabolic and efflux pump function;
- pharmacodynamics: the potential efficacy of a drug and whether the patient is predisposed to certain atypical effects.

These results do not change with age, but their interpretation can evolve as new data becomes available. Therefore, the Biron PGx reports are updated periodically. These results can also be useful for other medications, not covered by the report.

How to use pharmacogenomic recommendations

- 1. Only medications relevant for your patient need to be consulted.
- 2. Use the **Exposure** column to adjust doses for adequate plasma concentrations.
- 3. Use the Efficacy and Risk of atypical effect columns to choose the most compatible medication.

The **Exposure**, **Efficacy** and **Risk of atypical effect** columns are interpreted independently from each other. Medications are ordered by class with the most compatible options listed first within each class.

Efficacy


A higher dose may be required to achieve adequate When choosing between multiple clinically appropriate medications, you may give plasma concentrations. preference to a medication in which a lower number of variants have been identified (e.g., 2/2 is better than 4/6), in terms of their association with an increased likelihood of a A lower dose may be required to achieve adequate plasma poorer response or an atypical effect. concentrations. Signifies the presence of a high-impact Signifies the presence of variants Several metabolic pathways are involved, but their gene variant, which increases the associated with an increased risk of particular side effects, compared to capacities are opposed (e.g., PM and UM). Thus, a probability of a poorer response. calculation of dose adjustments is not possible based on non-carriers. current data and closer monitoring is recommended. Signifies that all of the tested variants predict an increased likelihood of a better X Medication not recommended by Drug not recommended by peer-reviewed guidelines due response, compared to non-carriers. This peer-reviewed guidelines due to a to a risk of toxicity or lack of efficacy. risk of severe side effects. medication may be a good option.

The notification "Genetic influence not available" signifies that there is currently no available data allowing for a genetically-based prediction of medication effect.

Nomenclature for enzyme phenotypes

(e.g., cytochrome P450s or CYP)

Exposure



NM is generally used to establish standard doses. This dose may be too high for **PM/IM** or too low for **RM/UM**, warranting a dose adjustments or the consideration of an alternative agent. For a pro-drug (e.g., clopidogrel, tramadol), phenotype variability will have the opposite effect.

<u>Inducible Metaboliser</u> (Ind) - Specific for CYP1A2, which can have increased function in the presence of an inducer, such as tobacco smoke, comparable to RM/UM.

Phone: 1-855-943-6379 Email: genetics@biron.com biron.com/en/genetics/pharmacogenomics



Risk of atypical effect

PHARMACOGENOMIC REPORT



Pain management

CONFIDENTIAL

To download the latest version, go here: secur.biron.com/login.

YOUR RESULTS ARE CONFIDENTIAL. As per the Genetic Non-Discrimination Act (S-201), no person, company or institution, including insurers and employers, can force you to share this report.

DO NOT MAKE ANY CHANGES TO YOUR CURRENT MEDICATION(S) WITHOUT TALKING TO YOUR DOCTOR FIRST. While genetics is important, other factors also contribute to how you react to medications. The final choice of medication used will be based on your health care provider's professional judgement and may be different than what is recommended in this report. This test does not determine your risk of any health problem. It only evaluates select portions of your DNA that help predict how you may react to the medications covered. For more information, visit biron.com/pgxtest.

	MINISTRATIVE D	ΑΙ	А
--	----------------	----	---

Patient Name Ordering Clinician Sample ID: BIO2107270327

Test-Firstname Test-Lastname Meredith Grey Sample Type: test

Patient Address Gender: Female Date ordered: 2022-08-21

1212 some street Date of birth: 1999-01-01 Date of sample reception: 2025-01-14

Ste-foy, Québec Phone Number: (418) 999-9999 Date of report: 2025-01-14 G2J4M5

Email: testsample.BIO2107270327@biron.local

Clinical Support

Email: genetique@biron.com Phone: 1-855-943-6379 Fax: (514) 317-2241

VARIANT(S) OF VERY IMPORTANT PHARMACOGENES (VIPS)

CYP1A2 IND, CYP2C9 IM, CYP2C19 IM, POR IA, CES1 IM, UGT1A1 IM, UGT1A4 IM, UGT2B7 variable, UGT2B15 IM.

Very Important Pharmacogenes (VIPs) - visit pharmgkb.org/vips for more information. Abbreviations - NM: Normal Metaboliser, IM: Intermediate Metaboliser, PM: Poor Metaboliser, RM: Rapid Metaboliser, UM: Ultrarapid Metaboliser, Ind: Inducible Metaboliser, NA: Normal Activity, IA: Intermediate Activity, PA: Poor Activity.

CAUTIONARY INFORMATION - MEDICATIONS TO AVOID OR USE WITH CAUTION

Medication	Identified risk	Recommendation
Carbamazepine Tegretol®	Risk of hypersensitivity reaction due to the presence of the HLA-B*15:02 allele.	If patient is carbamazepine-naïve, avoid carbamazepine. If the patient has previously used carbamazepine consistently for longer than three months without incidence of cutaneous adverse reactions, cautiously consider using carbamazepine. ¹
Oxcarbazepine Trileptal®	Risk of hypersensitivity reaction due to the presence of the HLA-B*15:02 allele.	If patient is oxcarbazepine-naïve, avoid oxcarbazepine. If the patient has previously used oxcarbazepine consistently for longer than three months without incidence of cutaneous adverse reactions, cautiously consider use of oxcarbazepine. ¹
Phenytoin Dilantin®	Risk of hypersensitivity reaction due to the presence of the HLA-B*15:02 allele.	If patient is phenytoin-naïve, avoid phenytoin. If the patient has previously used phenytoin consistently for longer than three months without incidence of cutaneous adverse reactions, cautiously consider use of phenytoïn.
Clopidogrel Plavix®	Reduced clopidogrel active metabolite formation; increased on-treatment platelet reactivity; increased risk for adverse cardiac and cerebrovascular events (CYP2C19 IM).	Cardiovascular indications: avoid standard dose (75mg/day) if possible; use prasugrel or ticagrelor at standard dose if no contraindication. Neurovascular:

CAUTIONARY INFORMATION - MEDICATIONS TO AVOID OR USE WITH CAUTION

Medication	ldentified risk	Recommendation
		consider alternative P2Y12 inhibitor at standard dose if clinically indicated and no contraindication. ²

PGx RECOMMENDATIONS - PSYCHIATRY

Dose increase may be required.

CONFIDENTIAL

Increased probability of a better response.

Dose reduction may be required.

- Greater potential for a poorer response or atypical effect.
- Exposure is difficult to predict, insufficient data to calculate dose
- Medication not recommended by peer-reviewed guidelines.

	G	enetic Associations Identific	ad .
Medications	Exposure	Efficacy	Risk of atypical effect
Antidepressants			
Selective serotoni	in reuptake inhibitors (SSRIs)		
Fluoxetine (Prozac®)	CYP2D6 NM, CYP2C9 IM - Initiate with recommended dose; a low dose may be adequate.	4/6 variants: increased likelihood of a poorer response (SLC6A4 - association demonstrated mainly in Caucasians, not observed in Asians, BDNF, FKBP5, HTR7); consider adding L-methylfolate 15mg/day (MTHFR).	Genetic influence not available
Fluvoxamine (Luvox®)	CYP2D6 NM, ABCB1 - Initiate with recommended dose but may require a higher dose.	4/6 variants: increased likelihood of a poorer response (SLC6A4 - association demonstrated mainly in Caucasians, not observed in Asians, BDNF, FKBP5, HTR7); consider adding L-methylfolate 15mg/day (MTHFR).	Genetic influence not available
Paroxetine (Paxil®)	CYP2D6 NM, ABCB1 - Initiate therapy with recommended starting dose but may require a higher dose.	4/6 variants: increased likelihood of a poorer response (SLC6A4 - association demonstrated mainly in Caucasians, not observed in Asians, BDNF, FKBP5, HTR7); consider adding L-methylfolate 15mg/day (MTHFR).	Genetic influence not available
Citalopram (Celexa®)	CYP2C19 IM - Initiate with recommended dose but consider a slower titration schedule and do not exceed the following daily doses: 30mg for adults up to 65 yrs; 15mg for adults 65 yrs or older. ^{3,4}	4/6 variants: increased likelihood of a poorer response (SLC6A4 - association demonstrated mainly in Caucasians, not observed in Asians, BDNF, FKBP5, HTR7); consider adding L-methylfolate 15mg/day (MTHFR).	1/1 variant: increased risk of gastrointestinal side effects (HTR2A).
Escitalopram (Cipralex®)	CYP2C19 IM - Initiate with recommended dose but consider a slower titration schedule and do not exceed the following daily doses: 15mg for adults up to 65 yrs; 7.5mg for adults 65 yrs or older. ^{3,4}	4/6 variants: increased likelihood of a poorer response (SLC6A4 - association demonstrated mainly in Caucasians, not observed in Asians, BDNF, FKBP5, HTR7); consider adding L-methylfolate 15mg/day (MTHFR).	1/1 variant: increased risk of gastrointestinal side effects (HTR2A).
Sertraline (Zoloft®)	CYP2B6 NM, CYP2C19 IM - Initiate therapy with recommended starting dose but with a slower titration schedule and a lower maintenance dose. ³	4/5 variants: increased likelihood of a poorer response (SLC6A4 - association demonstrated mainly in Caucasians, not observed in Asians, BDNF, FKBP5, HTR7); consider adding L-methylfolate 15mg/day (MTHFR).	1/1 variant: increased risk of gastrointestinal side effects (HTR2A).
Serotonin-norepi	inephrine reuptake inhibitors (SNRIs)		
Desvenlafaxine (Pristiq®)	UGT1A1 IM, UGT2B15 IM - Initiate with recommended dose; a low dose may be adequate.	1/2 variants: increased likelihood of a poorer response (FKBP5); consider adding L-methylfolate 15mg/day (MTHFR).	Genetic influence not available
Venlafaxine-XR (Effexor XR®)	CYP2D6 NM, ABCB1 - Initiate with recommended dose but may require a higher dose.	3/4 variants: increased likelihood of a poorer response (COMT, FKBP5, SLC6A2); consider adding L-methylfolate 15mg/day (MTHFR).	Genetic influence not available

PGx RECOMMENDATIONS - PAIN MANAGEMENT

Dose increase may be required.

CONFIDENTIAL

Increased probability of a better response.

Dose reduction may be required.

Greater potential for a poorer response or atypical effect.

Exposure is difficult to predict, insufficient data to calculate dose

Medication not recommended by peer-reviewed guidelines.

	Genetic Associations Identified						
Medications	Exposure	Efficacy	Risk of atypical effect				
Analgesic							
Acetaminophen (Tylenol®)	Genetic influence not available	Genetic influence not available	Genetic influence not available				
Antimetabolite							
Methotrexate	ABCB1, MTHFR - Consider using a lower dose.	Genetic influence not available	Genetic influence not available				
Cannabinoids							
Cannabidiol (CBD)	CYP3A4 NM, CYP2C9 IM - Normal exposure	Genetic influence not available	Genetic influence not available				
Nabilone (Cesamet®)	Genetic influence not available	Genetic influence not available	Genetic influence not available				
Tetrahydrocannabinol (THC)	CYP3A4 NM, CYP2C9 IM - Normal exposure	Genetic influence not available	1/2 variants: increased risk of cannabis use disorder (CNR1).				
Muscle Relaxant							
Carisoprodol (Soma®)	CYP2C19 IM - Initiate therapy with recommended starting dose and use with caution; a low dose may be adequate.	Genetic influence not available	Genetic influence not available				
Methocarbamol (Robaxin®)	Genetic influence not available	Genetic influence not available	Genetic influence not available				
Cyclobenzaprine (Flexeril®)	CYP3A4 NM, CYP1A2 Ind - Initiate with recommended dose but may require a higher dose, especially with CYP1A2 inducers, such as smoke.	Genetic influence not available	Genetic influence not available				
Tizanidine (Zanaflex®)	CYP1A2 Ind - Initiate with recommended dose but may require a higher dose, especially with CYP1A2 inducers, such as smoke.	Genetic influence not available	Genetic influence not available				
Nonsteroidal Anti-I	nflammatory Drugs (NSAID)						
Acetylsalicylic acid (Aspirin®)	Genetic influence not available	Genetic influence not available	Genetic influence not available				
Diclofenac (Voltaren®)	UGT2B7 IM - Initiate with recommended dose; a low dose may be adequate.	Genetic influence not available	Genetic influence not available				
Etodolac (Ultradol®)	CYP2C9 IM - Initiate with recommended dose; a low dose may be adequate.	Genetic influence not available	Genetic influence not available				
Indomethacin (Indocid®)	CYP2C9 IM - Initiate with recommended dose; a low dose may be adequate.	Genetic influence not available	Genetic influence not available				
Ketorolac (Toradol®)	CYP2C9 IM - Initiate with recommended dose; a low dose may be adequate.	Genetic influence not available	Genetic influence not available				
Naproxen (Naprosyn®)	Genetic influence not available	Genetic influence not available	Genetic influence not available				

	Genetic Associations Identified						
Medications	Exposure	Efficacy	Risk of atypical effect				
Hydromorphone (Dilaudid®)	Genetic influence not available	0/1 variant: no increased likelihood of a poorer response.	Genetic influence not available				
Meperidine (Demerol®)	CYP2B6 NM, POR, CYP3A4 NM - Initiate therapy with recommended starting dose; a low dose may be adequate.	Genetic influence not available	Genetic influence not available				
Methadone	CYP2B6 NM, CYP3A4 NM, POR - Initiate therapy with recommended starting dose; a low dose may be adequate.	Genetic influence not available	Genetic influence not available				
Nalbuphine (Nubain®)	Genetic influence not available	Genetic influence not available	Genetic influence not available				
Oxycodone (Supeudol®)	CYP3A4 NM, CYP2D6 NM - Normal exposure	0/1 variant: no increased likelihood of a poorer response.	Genetic influence not available				
Remifentanil (Ultiva®)	Genetic influence not available	Genetic influence not available	Genetic influence not available				
Sufentanil (Sufenta®)	CYP3A4 NM - Normal exposure	Genetic influence not available	Genetic influence not available				
Tapentadol (Nucynta®)	Genetic influence not available	Genetic influence not available	Genetic influence not available				
Tramadol (Ultram®)	CYP2D6 NM - Normal exposure	0/1 variant: no increased likelihood of a poorer response.	Genetic influence not available				
Codeine	CYP2D6 NM, UGT2B7 RM - Initiate with recommended dose but may require a higher dose.	0/1 variant: no increased likelihood of a poorer response.	Genetic influence not available				
Morphine	UGT2B7 RM - Initiate with recommended dose but may require a higher dose.	0/1 variant: no increased likelihood of a poorer response.	0/1 variant: no increased risk of gastrointestinal side effects (FAAH).				
Opioid antagonist	ts						
Naltrexone (Revia®)	Genetic influence not available	1/1 variant: increased likelihood of a poorer response when used in combination with bupropion for weight loss (ANKK1).	Genetic influence not available				
Naloxone (Narcan®)	Genetic influence not available	1/1 variant: increased likelihood of a poorer response (OPRM1).	Genetic influence not available				

PGx RECOMMENDATIONS - COMPLEMENTARY TREATMENTS

Dose increase may be required.

CONFIDENTIAL

Increased probability of a better response.

Dose reduction may be required.


Greater potential for a poorer response or atypical effect.

Exposure is difficult to predict, insufficient data to calculate dose adjustments.

Medication not recommended by peer-reviewed guidelines.

	G	Genetic Associations Identified			
Medications	Exposure	Efficacy	Risk of atypical effect		
Antiemetics					
Dimenhydrinate (Gravol®)	Genetic influence not available	Genetic influence not available	Genetic influence not available		
Granisetron (Kytril®)	CYP3A4 NM, CYP3A5 PM - Normal exposure	Genetic influence not available	Genetic influence not available		
Ondansetron (Zofran®)	CYP2D6 NM - Normal exposure	Genetic influence not available	Genetic influence not available		
Palonosetron (Aloxi®)	CYP2D6 NM - Normal exposure	Genetic influence not available	Genetic influence not available		
Proton pump inhib	itors (PPI)				
Esomeprazole (Nexium®)	Genetic influence not available	Genetic influence not available	Genetic influence not available		
Dexiansoprazole (Dexilant®)	CYP2C19 IM - Initiate therapy with standard dose but for chronic therapy (> 12 weeks) and efficacy achieved, consider reducing the daily dose by 50% and monitor for continued efficacy. ⁸	Genetic influence not available	Genetic influence not available		
Lansoprazole (Prevacid®)	CYP2C19 IM - Initiate therapy with standard dose but for chronic therapy (> 12 weeks) and efficacy achieved, consider reducing the daily dose by 50% and monitor for continued efficacy. ⁸	Genetic influence not available	Genetic influence not available		
Omeprazole (Losec®)	CYP2C19 IM - Initiate therapy with standard dose but for chronic therapy (> 12 weeks) and efficacy achieved, consider reducing the daily dose by 50% and monitor for continued efficacy. ⁸	Genetic influence not available	Genetic influence not available		
Pantoprazole (Pantoloc®)	CYP2C19 IM - Initiate therapy with standard dose but for chronic therapy (> 12 weeks) and efficacy achieved, consider reducing the daily dose by 50% and monitor for continued efficacy. ⁸	Genetic influence not available	Genetic influence not available		

METABOLISM	(D O D O O O O O O O O O O O O O O O O O		e enzymes listed below influence your exposure to active metabolite(s).	
	CYP450	POR	CES	UGT
Variants detected / tested	3/8	1/1	1/1	4/4



PGx ASSOCIATIONS - EXPOSURE

DISTRIBUTION	The level of activity of the efflux pumps listed below influence your exposure to medications at the site of action (e.g., central nervous system).	
	ABCB1	ABCG2
Variant(s) detected / tested	1/1	1/1

Efflux Pump		Phenotype		Genotype	Details
	Poor activity	Intermediate activity	Normal activity		
EFFLUX PUN	MPS (P-GLY	COPROTEIN, MU	JLT IDRUG RI	ESISTANCE PF	ROTEIN)
ABCB1	+	-		AAT/ACT	Poor activity (PA) Reduced efflux pump activity may result in higher substrate concentrations, increasing the likelihood of response but also increasing the risk of adverse drug reactions.
ABCG2		+		G/T	Intermediate activity (IA) Reduced efflux pump activity may result in higher substrate concentrations, increasing the likelihood of response but also increasing the risk of adverse drug reactions.

The variants listed below have been associated with a higher likelihood of a poorer response to medications, compared to non-carriers.

Poorer response more likely

Variant(s) detected / tested

CONFIDENTIAL

10/18

Poorer response more likely			
Gene	Reference SNP cluster ID	Variant detected	Associated medication(s)
ANKK1	rs1800497	yes	naltrexone
BDNF	rs6265	yes	citalopram, escitalopram, sertraline, fluoxetine, fluvoxamine, paroxetine
COMT	rs4680	yes	venlafaxine (AA)
FKBP5	rs4713916	yes	citalopram, desvenlafaxine, duloxetine, escitalopram, sertraline, fluoxetine, fluoxamine, mirtazapine, paroxetine, venlafaxine
GNB3	rs5443	yes	clonidine
HTR7	rs7905446	yes	citalopram, escitalopram, sertraline, fluoxetine, fluvoxamine, paroxetine
OPRM1	rs1799971	yes	naloxone
SLC6A2	rs2242446	yes	venlafaxine
SLC6A4	5-HTTLPR	yes	citalopram, escitalopram, sertraline, fluoxetine, fluvoxamine, paroxetine
TPH2	rs1487278	yes	amitriptyline, clomipramine, desipramine, doxepin, imipramine, mirtazapine, nortriptyline, trimipramine
ANKK1	rs1800497	no	bupropion
BDNF	rs6265	no	esketamine, ketamine
DRD3	rs963468	no	duloxetine
GRIK1	rs2832407	no	topiramate
GRIK4	rs1954787	no	citalopram, desvenlafaxine, duloxetine, escitalopram, fluoxetine, fluvoxamine, paroxetine, sertraline, venlafaxine
HTR2A	rs2770296	no	bupropion
HTR2A	rs6311	no	citalopram, escitalopram, fluoxetine, fluvoxamine, paroxetine
OPRM1	rs1799971	no	codeine, hydrocodone, hydromorphone, morphine, tramadol

PGx ASSOCIATIONS - ADVERSE DRUG REACTIONS

The variants listed below have been associated with a higher likelihood of specific adverse drug reactions, compared to non-carriers.

Increased risk of adverse drug reactions

Variant(s) detected / tested

CONFIDENTIAL

4/7

Increased risk of adverse drug reactions					
Gene	Reference SNP cluster ID	Variant detected	Associated medication(s)		
ANKK1	rs1800497	yes	weight gain: valproic acid		
CNR1	rs806380	yes	cannabis use disorder: THC		
HLA-B	*15:02	yes	cutaneous adverse reactions: carbamazepine, oxcarbazepine, phenytoin		
HTR2A	rs6311	yes	gastrointestinal side effects (nausea, vomiting): citalopram, escitalopram, sertraline		
FAAH	rs324420	no	gastrointestinal side effects (nausea, vomiting): morphine (A allele); cannabis use disorder: THC (AA)		
HLA-A	*31:01	no	cutaneous adverse reactions: carbamazepine		
SLC6A2	rs28386840	no	hypertension: esketamine		

ANALYTICAL RESULTS

The following analytical results were used to generate the pharmacogenomic interpretations found in this report. Technical limitations inherent with the methods used to produce these results may hinder the attribution of a definitive phenotype (see "TEST METHODOLOGY AND LIMITATIONS").

Genes	Variant Detail	ls (GRCH38.p12)	Result
ABCB1	rs1045642	chr7:87509329	A A
	rs2032582	chr7:87531302	AJC
	rs2032583	chr7:87531245	T T
ABCG2	rs2231142	chr4:88131171	G T
ADRA2A	rs1800544	chr10:111076745	C G
ANKK1	rs1800497	chr11:113400106	G G
BDNF	rs6265	chr11:27658369	C C
CACNG2	rs2283967	chr22:36567486	C C
CES1	rs71647871	chr16:55823658	CIT
CNR1	rs806380	chr6:88154934	AA
COMT	rs4680	chr22:19963748	AIG
CYP1A2	rs762551	chr15:74749576	CIA
	rs2069514	chr15:74745879	GİG
CYP2A6	rs1801272	chr19:40848628	AlA
	rs5031017	chr19:40843845	G G
	rs28399433	chr19:40850474	A A
CYP2B6	rs2279343	chr19:41009358	G G
	rs3745274	chr19:41006936	GIG
	rs28399499	chr19:41012316	T T
CYP2C9	rs1057910	chr10:94981296	CIA
	rs1799853	chr10:94942290	C C
	rs7900194	chr10:94942309 chr10:94949282-94949283	GIG
	rs9332131 rs9332239	chr10:94949282-94949285	A A C C
	rs28371685	chr10:94981224	CIC
	rs28371686	chr10:94981301	CIC
	rs72558187	chr10:94941958	ΤİΤ
	rs72558190	chr10:94947782	C C
CYP2C19	rs4244285	chr10:94781859	G A
	rs4986893	chr10:94780653	GIG
	rs6413438	chr10:94781858	CIC
	rs12248560	chr10:94761900	C C
	rs12769205 rs17884712	chr10:94775367 chr10:94775489	G A G G
	rs28399504	chr10:94762706	A A
	rs41291556	chr10:94775416	TIT
	rs56337013	chr10:94852738	clc
	rs72552267	chr10:94775453	GG
	rs72558186	chr10:94781999	ΤİΤ
CYP2D6	rs16947	chr22:42127941	G G
	rs1065852	chr22:42130692	GIG
	rs1135840	chr22:42126611	C C
	rs3892097	chr22:42128945	CIC
	rs5030655	chr22:42129084	A A
	rs5030656 rs5030862	chr22:42128174-42128178 chr22:42130668	A A C C
	rs5030865	chr22:42129033	CIC
	rs5030867	chr22:42127856	T T
	rs28371725	chr22:42127803	clc
	rs28371706	chr22:42129770	GG
	rs35742686	chr22:42128242	ΤİΤ
	rs59421388	chr22:42127608	clc
	rs72549354	chr22:42128815-42128817	DID
	rs774671100	chr22:42130555-42130755	GIG
	rs201377835 Gene Deletion	chr22:42129910	C C
	Gene Deletion Gene Duplication	n/a n/a	Not Detecte Not Detecte
	GGTE Dupitedtion	11/0	NOL DELECTE

Genes	Variant Deta	ils (GRCH38.p12)	Result
CYP3A4	rs4986907	chr7:99769804	CIC
	rs35599367	chr7:99768693	GIG
	rs55785340	chr7:99768360	A A
	rs67666821	chr7:99758184-99758188	DD
	rs72552799	chr7:99770165	C C
CYP3A5	rs776746	chr7:99672916	CIC
	rs10264272 rs28365083	chr7:99665212 chr7:99652613	C C G G
	rs28383479	chr7:99660516	CIC
	rs41303343	chr7:99652771	DID
	rs55817950	chr7:99676198	GIG
	rs56411402	chr7:99665237	TIT
DRD2	rs6275	chr11:113412755	G G
DRD3	rs963468	chr3:114144040	AIG
FAAH	rs324420	chr1:46405089	C C
FKBP5	rs4713916	chr6:35702206	AIG
GNB3	rs5443	chr12:6845711	C T
GRIK1	rs2832407	chr21:29595188	CIC
GRIK4	rs1954787	chr11:120792654	CIC
HLA-	rs1061235	chr6:29945521	AlA
A*31:01			
HLA- B*15:02	rs144012689	chr6:31355003	A A
HTR2A	rs6311	chr13:46897343	C C
	rs6313	chr13:46895805	GIG
	rs2770296	chr13:46866425	C C
HTR2C	rs3813929	chrX:114584047	C C
HTR7	rs7905446	chr10:90859404	G T
INSIG2	rs17047764	chr2:118111006	G G
long non-	rs74795342	chr21:18954018	AIG
coding (Inc) RNA	rs75222709	chr21:18955109	G T
MC4R	rs489693	chr18:60215554	AIC
	rs17782313	chr18:60183864	T T
MTHFR	rs1801131	chr1:11794419	GIG
000144	rs1801133	chr1:11796321	GIG
OPRM1 POR	rs1799971	chr6:154039662	AlA
SLC6A2	rs2868177	chr7:75960585	AIG
SLC6A2	rs5569 rs2242446	chr16:55697923 chr16:55656513	G G T T
	rs28386840	chr16:55652906	AA
SLC6A4	5-HTTLPR	chr17:30190154-30240133	SIS
SLC6A5	rs2298826	chr11:20638211	GIG
TH	rs2070762	chr11:2165105	AIG
TPH2	rs1487278	chr12:72007071	CIT
UGT1A1	rs4148323	chr2:233760498	GIG
	rs34815109	chr2:233760234-233760248	6 7
UGT1A4	rs2011425	chr2:233718962	G T
UGT2B7	rs7439366	chr4:69098620	T C
UGT2B15	rs1902023	chr4:68670366	C A

TEST METHODOLOGY AND LIMITATIONS

The Biron pharmacogenomic test for psychiatry and pain management is a MALDI-TOF-based single nucleotide primer extension genotyping test; laboratory developed and validated test (LDT), not approved by Health Canada. Nucleic acid amplification techniques may be subject to general interference by factors such as reaction inhibitors and low quality or quantity of extracted DNA. Factors influencing the amount and quality of extracted DNA include but are not limited to patient oral hygiene, collection technique and presence of dietary or microbial source of nucleic acids and nuclease. When present, these interferents typically yield no result rather than an inaccurate one. Risk of suboptimal DNA quantity or quality is significantly reduced by automated DNA extraction which uses chemistry without PCR inhibitors (magnetic beads) and systematic dilution, quantitation and normalization of DNA before nucleic acid amplification. Very infrequent variants or polymorphisms occurring in primer-binding regions may also affect testing and could produce an erroneous result or assay failure. The test does not detect all known and unknown variations in the genes tested, nor does absence of a detectable variant (typically reported as *1 for metabolic enzymes) rule out the presence of other, non-detected variants. The test detects CYP2D6 deletion and duplication but cannot differentiate duplication in the presence of deletion. CYP2D6 deletion and duplication assays can translate into equivocal phenotype results where a range of enzyme activity level must be reported. Test results and clinical interpretation may be inaccurate for individuals who have undergone or are receiving non-autologous blood transfusions, tissue, and/or organ transplant therapies.

DISCLAIMER

CONFIDENTIAL

Biron Health Group developed this pharmacogenomic report. This test does not diagnose any disorder, condition or disease. The interpretations and recommendations provided in this report are intended as a clinical support tool (DST) to be used solely by a healthcare professional. Treatment decisions for the patient remain the sole responsibility of the treating healthcare provider. The interpretations of the results provided by this report were determined by Biron's data curation protocol, which were established as per the current available scientific evidence available at the time this report version was created. As more evidence becomes available in the future, these interpretations may change. Some variants tested may not be used to provide report interpretations due to a lack of clear gene-drug association as determined by Biron's data curation protocol. The presence of a notification within the "Exposure", "Efficacy" or "Adverse Drug Reactions" categories for a given drug indicates that an associated genetic variant was detected. The lack of a notification within these categories for a given drug does not eliminate the requirement for dose adjustments for optimal dosage, does not guarantee effective drug therapy and does not eliminate the risks of adverse drug reactions. Commercial names are indicated as examples and do not consist an exhaustive list.

REFERENCES

For the full list of references, contact pgxinfo@biron.com

Reference(s) cited in this report:

- 1. Phillips EJ, et al. Clinical Pharmacogenetics Implementation Consortium Guideline for HLA Genotype and Use of Carbamazepine and Oxcarbazepine: 2017 Update. Clin Pharmacol Ther (2017).
- 2. Lee CR, et al. Clinical Pharmacogenetics Implementation Consortium Guideline for CYP2C19 Genotype and Clopidogrel Therapy: 2022 Update. Clin Pharmacol Ther (2022).
- 3. Bousman CA, et al. Clinical Pharmacogenetics Implementation Consortium (CPIC) Guideline for CYP2D6, CYP2C19, CYP2B6, SLC6A4, and HTR2A Genotypes and Serotonin Reuptake Inhibitor Antidepressants. Clin Pharmacol Ther (2023).
- 4. Brouwer J, et al. Dutch Pharmacogenetics Working Group (DPWG) guideline for the gene-drug interaction between CYP2C19 and CYP2D6 and SSRIs. Eur J Hum Genet (2021).
- 5. Dutch Pharmacogenetics Working Group (DPWG) 2021 update (https://api.pharmgkb.org/v1/download/file/attachment/DPWG_May_2021.pdf)
- Karnes JH, et al. Clinical Pharmacogenetics Implementation Consortium (CPIC) Guideline for CYP2C9 and HLA-B Genotypes and Phenytoin Dosing: 2020 Update. Clin Pharmacol Ther (2021).
- 7. Theken KN, et al. Clinical Pharmacogenetics Implementation Consortium Guideline (CPIC) for CYP2C9 and Nonsteroidal Anti-Inflammatory Drugs. Clin Pharmacol Ther (2020).
- 8. Lima JJ, et al. Clinical Pharmacogenetics Implementation Consortium (CPIC) Guideline for CYP2C19 and Proton Pump Inhibitor Dosing. Clin Pharmacol Ther (2021).