



THE ONLY POST-VASECTOMY SPERMOGRAM USING FLOW CYTOMETRY

Eliminates chance
and the subjectivity
of microscopic
observation



The best test on the market
for confirming the
absence of sperm
in a specimen



NEW NASCI POST- VASECTOMY SPERMOGRAM

**Much more reliable
and convenient**
for patients



To learn more, visit biron.ca


THE MOST ACCURATE SPERM COUNT IN CANADA

A sperm count below 5,000 sperm/ml is generally indicative of a successful vasectomy and no repeat test should be necessary.

Analysis of a sperm volume that is
100 to 1,000x
larger

Using a technique that is
10x more sensitive

For more than 80% of patients,
the success of the vasectomy
can be confirmed right
from the first spermogram.



In the past, specimens had to be kept at body temperature and analyzed within 2 hours of collection:

30%
of men

did not do the test
(instructions too
inconvenient or laboratory
too far away).

There was a lot of insecurity
and anxiety about the accuracy
of the results.

To learn more,
visit biron.ca


Biron
Laboratoire médical



BENEFIT FROM THE RELIABILITY OF THE NASCI POST-VASECTOMY SPERMOGRAM

Opt for a test that is:

- ✓ more reliable
- ✓ more convenient
- ✓ available closer to home